

## Which type of Masks to choose?

ASTM Level 1 Mask	ASTM Level 2 Mask	ASTM Level 3 Mask	EN Type II Mask	N95
ASTM2100-11	ASTM2100-11	ASTM2100-11	EN14638:2014	NOISH & FDA
Protection for Droplet Virus	Protection for Droplet Virus	Protection for Droplet Virus	Protection for Droplet Virus	Protection for Air Borne Virus
@80mm Hg ≥95% PFE @0.1 micron ≥95% of BFE @3 micron <4.0mm H2O/cm2 Class 1 Flame Spread	@120mm Hg ≥98% PFE @0.1 micron ≥98% of BFE @3 micron <5.0mm H2O/cm2 Class 1 Flame Spread	@160mm Hg ≥98% PFE @0.1 micron ≥98% of BFE @3 micron <5.0mm H2O/cm2 Class 1 Flame Spread	N/A N/A ≥98% of BFE @3 micron <1.5mm H2O/cm2 N/A	@160mm Hg 99.9% PFE 0.1 micron 95% of BFE 3 micron >5.0mm H2O/cm2 Class 1 Flame Spread
Procedure Mask 2011-2018M, 2120, 2331 2070, 2330, 2000B, 2150, 2026	Surgical Mask 2882, 2885, 2879, 2880 2127, 2010SBL, 2010SWH 2010SPK 2010SYL, KMI3531120	Surgical Mask KMI34301F1 2025, 2050BF	Procedure Mask 2394	N95 2321, 9500